

New approach of simple pretreatment  
for GCMS analysis  
by using monolithic material  
sorptive extraction

- **Current sample pretreatment techniques**
- **Monolithic material sorptive extraction (MMSE)**
- **PAHs analysis with MMSE in water sample**
- **Conclusion**



# Introduction

- Environmental pollutions and food contaminations are biggest social concerns, nowadays. Because of the complex matrices in environmental substances and food stuffs, sample preparations for the quantitative analysis of toxic and hazardous compounds with chromatographic instruments are inevitable in order to obtain reliable results. Gas chromatography mass spectrometry (GCMS) has become the important means in the areas of qualitative and quantitative analysis of environmental monitoring and trace amount of contaminants in food stuffs.
- Although the performances of sensitivity and selectivity of this instrument have rapidly been improved in recent years, the sample pretreatment procedures remain the bottleneck for rapid analysis. Therefore, economical and effective pretreatment methods to eliminate matrix interferences, to achieve high throughput and to enrich the analyte for high sensitive detection are demanded.

# Goal for sample pretreatment

There are many types of sample pretreatment techniques. However, most of the goal for sample pretreatment are the same.

- High extraction efficiency
- Removal of potential interferences
- Enrichment of an analyte concentration
- Derivatization of an analyte for better separation and higher sensitivity

Better accuracy and precision must be taken into account in sample pretreatment techniques.

## Current sample pretreatment techniques

- Soxhlet extraction
- Liquid-liquid extraction (LLE)
- Ultrasonic extraction
- Solid phase extraction (SPE)
- Solid phase microextraction (SPME)
- Stir-bar sorptive extraction (SBSE)
- **Monolithic material sorptive extraction (MMSE)**
- Matrix solid-phase dispersion (MSPD)
- Liquid phase microextraction (LPME)
- Supercritical fluid extraction (SFE)
- Pressurized liquid extraction (PLE)
- Microwave-assisted extraction (MASE)
- Superheated water extraction (SWE)

## Trends of sample pretreatment



Small initial sample sizes

Small preparation volumes

Use of less toxic organic solvents or no organic solvents

Introduction of automation to eliminate human error.


# Silica Gel Monolithic Materials



# Monolithic Material Sorptive Extraction (MMSE)

Because of through-pore, not only **gas substances** but also **liquid substances** are permeable. Meso-pore enlarges the surface area to contact substances. This material can be also used as a sample pretreatment media. This technique may be same as stir bar sorptive extraction (SBSE) and solid-phase microextraction (SPME). However, the difference between MMSE and SBSE/SPME is surface area which will be affected to **high extraction efficiency** resulting in **high sensitivity** (Surface area:150 m<sup>2</sup>/g).

# Type and Dimension of MMSE



## Type

(1) ODS (white)

(2) ODS with active carbon (black)

(3) ODS with graphite carbon (gray)

## Dimension

Disk type: 10 mmOD x 1 mm thickness

Rod type: 2.9 mmOD x 5 mm height



# Selection of MMSE

**Thermal Desorption**

**Solvent Extraction**




ODS with graphite carbon (gray)



ODS (white)




ODS with active carbon (black)



**Adsorption capability**

**Thermal Desorption capability**

# Target Application



Gases Sample

Liquid Sample

Fragrance & Flavor

Bad Smell

Environmental  
-PAHs  
-Pesticides

Plant Flower

Food, Drink, Wine

# Sampling with MMSE

Passive Gas Sampling



Head Space Gas sampling



Agitating Liquid sampling



# Agitation for liquid sample



Incubator



Shaker

# Desorption from MMSE

- Solvent extraction



Add  
organic  
solvent



Very low  
volume



Ultra-  
sonication  
for 5 min



- Thermal desorption




Thermal desorption device for direct injection

No organic solvents

High recovery, high sensitivity


# Thermal desorption device for MMSE



# Coffee flavor analysis with MMSE vs. SPME



(coffee)




Thermal desorption for MMSE

# Peach juice analysis with MMSE vs. SBSE

(SBSE:Stir-Bar Sorptive Extraction)



Peach juice




MMSE


SBSE

- #1;beta-Linalool
- #2;gamma-Caprolactone
- #3;alpha-Terpineol
- #4;Benzyl acetate

# Comparison of MMSE, SBSE and SPME - Log P and Recovery -



SAMPLE: 250ppb-aroma mix  
in 15% NaCl



temp: 60°C  
time: 30 min  
rpm: 90r.p.m

MonoTrap; MeCl<sub>2</sub> extraction




SPME; TD




SBSE; ACN; extraction

# Flavor analysis in candy



# Chromatogram of Candy flavor





- |   |                    |    |                       |
|---|--------------------|----|-----------------------|
| 1 | <i>β</i> -Myrcene  | 8  | Ethyl laevulinate     |
| 2 | Limonene           | 9  | <i>α</i> -Terpineol   |
| 3 | cis-3-Hexenol      | 10 | Citral                |
| 4 | Furfural           | 11 | <i>β</i> -Damascenone |
| 5 | Decanal            | 12 | Benzyl Alcohol        |
| 6 | <i>β</i> -Linalool | 13 | Phenylethyl Alcohol   |
| 7 | Propylene Glycol   |    |                       |

※Identified with spectral libraries

Red characters •••Identified with standard sample

# Flavor analysis in Japanese red wine



③ MT Extract Cup with Vial (Start UP-Kit)

## GC/MS analysis

**System** : SHIMADZU GC-2010, GCMS-QP2010

**Column** : **InertCap Pure-WAX**(Cat.1010-68142)**New!!**

0.25mm I.D. × 30m  $d_f=0.25 \mu\text{m}$

**Column Temp** : 40°C(5min) → 6°C/min → 250°C(5min)


**Carrier Gas** : He 95kPa

**Injection** : Split / Splitless, 1  $\mu\text{L}$

250°C

**Detection** : MS Scan (m/z; 55-400)

# Chromatogram of Japanese red wine flavor



# Identification of compounds in Japanese red wine

1	2,2,6-Trimethyl-6-vinyltetrahydropyran	18	Benzaldehyde
2	Isoamyl acetate	19	3-Ethyl-4-methylpentanol
3	Limonene	20	2-Bornene
4	1-Pentanol	21	n-Propyl propionate
5	Ethyl hexanoate	22	Ethyl dl-2-hydroxycaproate
6	Maleic anhydride	23	$\beta$ -Cyclocitral
7	3-Methylpentanol	24	Ethyl decanoate
8	1,1-Dimethoxy-2-propanol	25	$\alpha$ -D-Galactopyranose methyl glycoside
9	Ethyl 2-hexenoate	26	Diethyl succinate
10	1-Hexanol	27	3-(Methylthio)-1-propanol
11	cis-3-Hexen-1-ol	28	1,5,8-Trimethyl-1,2-dihydronaphthalene
12	Nonanal	29	Hexanoic acid
13	cis-2-Hexen-1-ol	30	Benzyl Alcohol
14	Ethyl 2-hydroxy-3-methylbutanoate	31	Phenylethyl Alcohol
15	Ethyl octanoate	32	Diethyl dl-malate
16	Furfural	33	Octanoic Acid
17	2-Ethyl-1-hexanol		

※Identified with spectral libraries

Red characters ... Referred to Japan Perfumery & Flavoring Association Encyclopedia of Flavor[food]

# Analytical conditions for PAHs analysis

Thermal desorption device : ATAS-GL OPTIC3 Linex

Desorption Temperature : 250 °C, Time 10min

Injection Temperature: 250 °C

Split ratio : 1:10

Carrier Gas : He (1 mL/min)

GC/MS instrument : Shimadzu GCMS-QP2010 Plus


Column : InertCap Pesticide 0.25 mmIDx30 m

Column Temperature : 50 °C (3 min) → 10 °C/min → 310 °C

Detection : SIM mode (I.F.; 340 °C, I.S.; 230 °C)

# Chromatogram of STD PAHs

Peak #	Name
1	Naphthalene
2	Acenaphthylene
3	Acenaphthene
4	Fluorene
5	Phenanthrene
6	Anthracene
7	Fluoranthene
8	Pyrene
9	Benz[a]anthracene
10	Chrysene
11	Benzo[b]fluoranthene
12	Benzo[k]fluoranthene
13	Benzo[a]pyrene
14	Indeno[1,2,3-cd]pyrene
15	Dibenz[a,h]anthracene
16	Benzo[g,h,i]perylene




# Evaluation of sampling time for PAHs analysis

Recovery

MonoTrap RSC18 x 2

(0.025ppb - 0.67ppb)



- ◆ Naphthalene
- Acenaphthylene
- ▲ Acenaphthene
- Fluorene
- ◆ Phenanthrene
- Anthracene
- Fluoranthene
- Pyrene
- Benz[a]anthracene
- ◆ Chrysene
- Benzo[b]fluoranthene
- ▲ Benzo[k]fluoranthene
- Benzo[a]pyrene
- ▲ Benzo[g,h,i]perylene
- Dibenz[a,h]anthracene
- Indenno[1,2,3-cd]pyrene

# Recovery result

	<b>1 Rod</b>	<b>2 Rods</b>
Naphthalene	30%	51%
Acenaphthylene	51%	72%
Acenaphthene	65%	81%
Fluorene	70%	90%
Phenanthrene	79%	105%
Anthracene	84%	100%
Fluoranthene	84%	113%
Pyrene	81%	109%
Benz[a]anthracene	86%	113%
Chrysene	87%	108%
Benzo[b]fluoranthene	88%	117%
Benzo[k]fluoranthene	86%	110%
Benzo[a]pyrene	82%	112%
Benzo[g,h,i]perylene	85%	104%
Dibenz[a,h]anthracene	87%	102%
Indeno[1,2,3-	73%	87%

Sampling time : 3 hr

# Repeatability (0.25 ppb)

PAHs	n1	n2	n3	Ave	RSD
Naphthalene	214862	249466	248469	237599.0	8.3%
Acenaphthylene	232528	257623	280506	256885.7	9.3%
Acenaphthene	181728	200323	215534	199195.0	8.5%
Fluorene	200912	215331	229094	215112.3	6.6%
Phenanthrene	290295	309267	318268	305943.3	4.7%
Anthracene-d10(IS)	21452	21926	23885	22421.0	5.8%
Anthracene	274284	266254	254799	265112.3	3.7%
Fluoranthene	287263	302767	307833	299287.7	3.6%
Pyrene	286815	298771	310061	298549.0	3.9%
Benz[a]anthracene	214072	207132	197872	206358.7	3.9%
Chrysene-d12(IS)	20330	19863	23610	21267.7	9.6%
Chrysene	213968	209389	190320	204559.0	6.1%
Benzo[b]fluoranthene	226864	211373	205630	214622.3	5.1%
Benzo[k]fluoranthene	226928	215212	192301	211480.3	8.3%
Benzo[a]pyrene	157758	148425	143156	149779.7	4.9%
Benzo[g,h,i]perylene	187547	176735	168432	177571.3	5.4%
Dibenz[a,h]anthracene	192653	187190	171664	183835.7	5.9%
Indenno[1,2,3-cd]pyrane	188182	186530	175394	183368.7	3.8%

# Linearity result


Concentration (ppb)	0.1	1	10	100	r <sup>2</sup>
Naphthalene	8810	85133	908151	5548551	0.9961
Acenaphthylene	7848	76464	872538	6703493	0.9991
Acenaphthene	10529	103589	1107188	7337325	0.9975
Fluorene	2776	25793	272269	2365694	0.9998
Phenanthrene	4257	32143	305547	2496179	0.9996
Anthracene	625	7815	66312	597123	0.9999
Fluoranthene	7315	79293	778395	5450872	0.9983
Pyrene	8362	92965	919141	6139653	0.9977
Benz[a]anthracene	2773	36124	403727	3024379	0.9989
Chrysene	2697	35329	366602	2605325	0.9984
Benzo[b]fluoranthene	2718	38055	416201	3132210	0.9990
Benzo[k]fluoranthene	3457	47304	494551	3577392	0.9986
Benzo[a]pyrene	2962	37944	441353	3505523	0.9993
Benzo[g,h,i]perylene	1682	31334	392890	2991132	0.9990
Dibenz[a,h]anthracene	1852	26736	324509	1520019	0.9872
Indeno[1,2,3-cd]pyrene	1757	30934	352939	2091641	0.9954

# Limit of Quantitation (LOQ)


PAHs	LOQ (ppb)	PAHs	LOQ (ppb)
Naphthalene	0.065	Benz[a]anthracene	0.013
Acenaphthylene	0.050	Chrysene	0.012
Acenaphthene	0.067	Benzo[b]fluoranthene	0.013
Fluorene	0.015	Benzo[k]fluoranthene	0.015
Phenanthrene	0.011	Benzo[a]pyrene	0.016
Anthracene	0.003	Benzo[g,h,i]perylene	0.014
Fluoranthene	0.025	Dibenz[a,h]anthracene	0.011
Pyrene	0.027	Indenno[1,2,3-cd]pyrane	0.012

(LOQ: S/N=10)


# Chromatogram of Anthracene



# Measurement of PAHs in river water by MMSE



# Analysis of the bad breath after eating a garlic



## GC/MS

System

: SHIMADZU GC-2010、GCMS-QP2010

Column

: InertCap Pure-WAX

0.25mmI.D. × 60m df=0.25µm

Column Temp

: 35°C(5min)→4°C/min→250°C(10min)

Carrier Gas : He 120kPa

Injection


: Split less、1µL 250°C

Detection


: MS Scan (m/z;35-500)



1, Allyl Methyl sulfide





2, Allyl sulfide




3, Diallyl disulfide

# Result of breath analysis



Allyl sulfide



Diallyl disulfide

The main aroma compound of a raw garlic are almost decomposed or metabolized in a body, and are considered that allyl methyl sulfide is generated.

## Breath analysis with MMSE



**Breath can indicate some of symptoms of diseases.  
Collecting the breath with MMSE is extremely simple and the thermal desorption  
GCMS technique can easily determine organic substances as a maker of health  
conditions.**

## Conclusions

- MMSE can be used for simple pretreatment of the sample as a qualitative and quantitative analysis for gases and liquid samples.
- Small quantity of organic solvent is required for organic solvent desorption from monolithic materials.
- No organic solvent is required for thermal desorption from monolithic materials.
- Thermal desorption GCMS combined with MMSE has been confirmed as a simple pretreatment method for determination of 16 PAHs with acceptable analytical performances.
- As compared to SBSE and SMPE, MMSE had better sensitivity because of higher recovery.