


# Extraction of Ethyl Glucuronide (EtG) from Urine Using EVOLUTE AX 96-Well Plates and Columns

## Introduction

This application note describes the extraction of ethyl glucuronide (EtG) from urine using EVOLUTE AX plates or columns.



**Figure 1.** Ethyl glucuronide

Ethyl glucuronide (EtG, **Figure 1.**) is an ethanol metabolite formed by the glucuronidation of the parent molecule. The presence of EtG in urine and plasma is a useful and sensitive biomarker of alcohol intake in forensic toxicology. EtG can be detected at very low levels and being a metabolite it is useful in distinguishing between ingested ethanol and that which is adsorbed through the skin from ethanol containing products such as alcohol based hand washes. Rapid and reliable methods with robust extraction protocols are essential for analysis and quantitation of EtG from various matrices.

EVOLUTE AX is a resin-based mixed-mode sorbent with an optimised combination of non-polar (hydrophobic), polar (hydrophilic) and strong anion exchange interactions ideally constructed to extract acidic compounds like ethyl glucuronide.

## Sample Preparation Procedure

### Analytes

Ethyl glucuronide (EtG)

### EVOLUTE AX Configuration

25 mg 96-fixed-well plate (part number 603-0025-P01)

This application note outlines the procedure using the 25mg 96-well plate format. Other formats are available and optimized methodology is detailed in Table 1.

### EVOLUTE AX Procedure

Sample:	Urine (50 $\mu$ L ) spiked at 500 ng/mL with EtG
Sample Pre-treatment:	Urine (50 $\mu$ L) in acetonitrile (950 $\mu$ L)
Column Conditioning:	Methanol (1 mL)
Column Equilibration 1:	Water (1 mL)
Column Equilibration 2:	Acetonitrile (1 mL)
Sample Loading:	Pre-treated sample (1 mL)
Interference Elution 1:	Water (1 mL)
Interference Elution 2:	Acetonitrile (1 mL)
Analyte Elution:	Acetonitrile:water:formic acid, (1 mL, 95:4:1, v/v)
Post Extraction:	Evaporate to dryness and reconstitute with water/methanol (200 $\mu$ L, 80:20, v/v) for LC-MS/MS analysis.

This method is suitable to be applied across the range of EVOLUTE AX formats, dilution factors have been optimized to maximize recoveries and minimize ion suppression.

<b>EVOLUTE AX procedure</b>	<b>Solution</b>	25 mg 96-Well-Plate <b>Part No.</b> 603-0025-P01	50 mg/3 mL Column <b>Part No.</b> 613-0005-B	100 mg/3 mL Column <b>Part No.</b> 613-0010-B	200 mg /3 mL Column <b>Part No.</b> 613-0020-B	500 mg/6 mL Column <b>Part No.</b> 613-0050-C
Sample	Urine	50 µL	200 µL	200 µL	400 µL	1 mL
Sample pre-treatment	Acetonitrile	950 µL	1.8 mL	1.8 mL	3.6 mL	4 mL
Column Conditioning	Methanol	1 mL	2 mL	3 mL	6 mL	6 mL
Column equilibration 1	Water	1 mL	2 mL	3 mL	6 mL	6 mL
Column equilibration 2	Acetonitrile	1 mL	2 mL	3 mL	6mL	6 mL
Sample loading	Pre-treated sample	1 mL	2 mL	2 mL	4 mL	5 mL
Interference elution 1	Water	1 mL	2 mL	3 mL	6 mL	6 mL
Interference elution 2	Acetonitrile	1 mL	2 mL	3 mL	6 mL	6 mL
Analyte elution	Acetonitrile:Water: Formic acid, (95:4:1, v/v)	1 mL	2 mL	3 mL	6 mL	6 mL
Reconstitution	Water/Methanol (80:20, v/v)	200 µL	200 µL	200 µL	200 µL	200 µL

**Table 1.** Optimized methodology for all EVOLUTE AX formats.

#### HPLC Conditions


Instrument: Waters Acquity UPLC (Waters Assoc., Milford, MA, USA).  
Column: Acquity UPLC BEH C18 column (1.7 µm, 50 x 2.1 mm id) (Waters Assoc., Milford, MA, USA).  
Mobile Phase: Isocratic, 0.1% (v/v) formic acid aq/methanol (80:20, v/v) at a flow rate of 0.5 mL/min.  
Injection Volume: 10 µL  
Column Temperature: 40 °C

#### Mass Spectrometry Conditions


Instrument: Premier XE triple quadrupole mass spectrometer (Waters Assoc., Manchester, UK) equipped with an electrospray interface for mass analysis. Negative ions were acquired in the multiple reaction monitoring mode. (MRM); 221.0 > 84.9 (quantifier), 221.0 > 74.9 (qualifier).  
Desolvation Temperature: 450 °C  
Ion Source Temperature: 150 °C  
Collision Gas Pressure: 3.46 x 10<sup>-3</sup> mbar

#### Results

All results show recoveries greater than 80% with RSD's below 10%. The method shows that optimization is scalable for sample volumes ranging from 50 µL up to 1 mL.



**Figure 2.** EtG extraction recoveries from urine using EVOLUTE AX



**Figure 3.** Chromatogram showing extracted EtG (500 ng/mL) using EVOLUTE AX 25 mg 96 well-plate

### References

This application note is based on the poster 'Extraction of Ethyl glucuronide (EtG) using a New Resin-Based Mixed-Mode Strong Anion Exchange SPE Sorbent Prior to LC-MS/MS Analysis', L Williams et al, presented at TIAFT, Bonn, August 29 - September 2, 2010.

### NORTH AMERICA

Main Office: +1 704 654 4900  
 Toll Free: +1 800 446 4752  
 Fax: +1 704 654 4917  
 Order Tel: +1 704 654 4900  
 press (4) at the auto attendant  
 Order Fax: +1 434 296 8217  
 ordermailbox@biotage.com  
 1-pointsupport@biotage.com

### EUROPE

Main Office: +46 18 56 5900  
 Fax: +46 18 59 1922  
 Order Tel: +46 18 56 57 10  
 Order Fax: +46 18 56 57 05  
 order@eu.biotage.com

### Japan

Tel: +81 3 5627 3123  
 Fax: +81 3 5627 3121  
 jp\_order@biotage.com

### Distributors

Please visit our Web site at  
[www.biotage.com](http://www.biotage.com)  
 for contact details.