

## How to Prepare Solutions and Buffers

### **1.0 M Acetic Acid:**

To 400 mL D.I. H<sub>2</sub>O add 28.6 mL glacial acetic acid. Dilute to 500 mL with D.I. H<sub>2</sub>O. Storage: 25 °C in glass or plastic. Stability: 6 months

### **100 mM Acetic Acid:**

Dilute 40 mL 1.0 M acetic acid to 400 mL with D.I. H<sub>2</sub>O. Mix. Storage: 25 °C in glass or plastic. Stability: 6 months

### **100 mM Acetate Buffer (pH 4.5):**

Dissolve 2.93 g sodium acetate trihydrate in 400 mL D.I. H<sub>2</sub>O; add 1.62 mL glacial acetic acid. Dilute to 500 mL with D.I. H<sub>2</sub>O. Mix. Adjust pH to 4.5 ± 0.1 with 100 mM sodium acetate or 100 mM acetic acid. Storage: 25 °C in glass or plastic. Stability: 6 months; Inspect daily for contamination.

### **1.0 M Acetate Buffer (pH 5.0):**

Dissolve 42.9 g sodium acetate trihydrate in 400 mL D.I. H<sub>2</sub>O; Add 10.4 mL glacial acetic acid. Dilute to 500 mL with D.I. H<sub>2</sub>O. Mix. Adjust pH to 5.0 ± 0.1 with 1.0 M sodium acetate or 1.0 M acetic acid. Storage: 25 °C in glass or plastic. Stability: 6 months; Inspect daily for contamination.

### **100 mM Acetate Buffer (pH 5.0):**

Dilute 40 mL 1.0 M acetate buffer to 400 mL with D.I. H<sub>2</sub>O. Mix. Storage: 25 °C in glass or plastic. Stability: 6 months

### **7.4 M Ammonium Hydroxide:**

To 50 mL D.I. H<sub>2</sub>O add 50 mL concentrated NH<sub>4</sub>OH. Mix. Storage: 25 °C in glass or fluoropolymer plastic. Stability: Storage condition dependent.

### **100 mM Hydrochloric Acid:**

To 400 mL D.I. H<sub>2</sub>O add 4.2 mL concentrated HCl. Dilute to 500 mL with D.I. H<sub>2</sub>O. Mix. Storage: 25 °C in glass or plastic. Stability: 6 months

### **Methanol /Ammonium Hydroxide (98:2):**

To 98 mL CH<sub>3</sub>OH add 2 mL concentrated NH<sub>4</sub>OH. Mix. Storage: 25°C in glass or fluoropolymer plastic. Stability: 1 day.

### **0.35 M Sodium Periodate:**

Add 37.5 g sodium periodate to a 500 mL volumetric flask, dilute to volume with D.I. H<sub>2</sub>O. Mix. Stability: 2 mos. at room temperature.

### **CH<sub>2</sub>Cl<sub>2</sub> / IPA / NH<sub>4</sub>OH (78:20:2):**

To 20 mL IPA, add 2 mL concentrated NH<sub>4</sub>OH. Mix. Add 78 mL CH<sub>2</sub>Cl<sub>2</sub>. Mix. Storage: 25 °C in glass or fluoropolymer plastic. Stability: 1 day

**100 mM Phosphate Buffer (pH 6.0):**

Dissolve 1.70 g  $\text{Na}_2\text{HPO}_4$  and 12.14 g  $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$  in 800 mL D.I.  $\text{H}_2\text{O}$ . Dilute to 1000 mL using D.I.  $\text{H}_2\text{O}$ . Mix. Adjust pH to  $6.0 \pm 0.1$  with 100 mM monobasic sodium phosphate (lowers pH) or 100 mM dibasic sodium phosphate (raises pH). Storage:  $5^\circ\text{C}$  in glass. Stability: 1 month; Inspect daily for contamination.

**500 mM Phosphoric Acid:**

To 400 mL D.I.  $\text{H}_2\text{O}$  add 17.0 mL concentrated phosphoric acid. Dilute to 500 mL with D.I.  $\text{H}_2\text{O}$ . Mix. Storage:  $25^\circ\text{C}$  in glass or plastic. Stability: 6 months

**1.0 M Sodium Acetate:**

Dissolve 13.6 g sodium acetate trihydrate in 90 mL D.I.  $\text{H}_2\text{O}$ . Dilute to 100 mL with D.I.  $\text{H}_2\text{O}$ . Mix. Storage:  $25^\circ\text{C}$  in glass or plastic. Stability: 6 months

**100 mM Sodium Acetate:**

Dilute 10 mL 1.0 M sodium acetate to 100 mL with D.I.  $\text{H}_2\text{O}$ . Mix. Storage:  $25^\circ\text{C}$  in glass or plastic. Stability: 6 months

**100 mM Sodium Borate:**

Dissolve 3.81 g  $\text{Na}_2\text{B}_4\text{O}_7 \cdot 10 \text{H}_2\text{O}$  in 90 mL D.I.  $\text{H}_2\text{O}$ . Dilute to 100 mL with D.I.  $\text{H}_2\text{O}$ . Mix. Storage:  $25^\circ\text{C}$  in glass or plastic. Stability: 6 months.

**100 mM Sodium Phosphate Dibasic:**

Dissolve 2.84 g  $\text{Na}_2\text{HPO}_4$  in 160 mL D.I.  $\text{H}_2\text{O}$ . Dilute to 200 mL using D.I.  $\text{H}_2\text{O}$ . Mix. Storage:  $5^\circ\text{C}$  in glass. Stability: 1 month; Inspect daily for contamination.

**100 mM Sodium Phosphate, Monobasic:**

Dissolve 2.76 g  $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$  in 160 mL D.I.  $\text{H}_2\text{O}$ . Dilute to 200 mL with D.I.  $\text{H}_2\text{O}$ . Mix. Storage:  $5^\circ\text{C}$  in glass. Stability: 1 month. Inspect daily for contamination.

**100 mM Sulfuric Acid:**

To 400 mL D.I.  $\text{H}_2\text{O}$  add 2.7 mL concentrated  $\text{H}_2\text{SO}_4$ . Dilute to 500 mL with D.I.  $\text{H}_2\text{O}$ . Mix. Storage:  $25^\circ\text{C}$  in glass or plastic. Stability: 6 months