



OVERVIEW / INTRODUCTION

Since the introduction and implementation of LC/MS as a staple analytical tool in forensic laboratories, there have been new approaches to sample preparation. The higher sensitivity of LC/MS and the ability to inject 'aqueous' containing samples directly into the instrument has opened new options for conventional sample preparations. The need for rapid turnaround time for a larger list of drugs has also put pressure on laboratories to find alternatives to traditional methods. The usual liquid-liquid and solid phase extraction processes have seen a growth of 'crash and shoot' or 'dilute and shoot' sample preparation methods. Although these latter methods work most of the time for certain applications (i.e. primarily urine samples), these alternatives have also introduced new shortcomings.

LC/MS analysis is very prone to matrix suppression phenomenon. The 'crash' or 'dilute' methods no longer remove matrix and concentrate samples but instead dilute the final eluate. These methods can raise the LOD and by definition, lower the sensitivity of the method. The diluted samples will still contain unwanted matrix that when introduced into the system can contaminate the instrumentation. In addition, these methods usually require a 10-15 minute centrifugation of the samples prior to injection. This step is done to eliminate any particulates that might get caught in either the guard column or more expensive LC columns. Most LC column packing particle sizes are not greater than 5µm and can therefore be subject to clogging by certain samples.

This poster describes a method that uses positive pressure and a solid phase sorbent bed built with small pore frits to quickly and efficiently prepare samples for LC/MS analysis. This method referred to as FAST, eliminates the timely centrifugation, reduces matrix suppression effects and removes particulates greater than ~ 1µm. Samples can be diluted at a ratio as low as 1:1, which is useful for analytes at very low concentrations.

METHOD

I. FAST Method – Opiates

Sample Dilution Ratio	Sample* Volume	Dilution** Volume
1:1	500 µL	500 µL
1:4	200 µL	800 µL
1:9	100 µL	900 µL

* If sample is hydrolyzed add appropriate aliquot volume after hydrolysis is complete.
** Diluent is 50:50 (Methanol: Distilled Water)

- Sample and diluents are added in an appropriately labeled tube. Add appropriate volume internal standard(s). It is recommended to use an internal standard volume of no more than 200 µL.
- Set up extraction manifold with FAST cartridges and auto-sampler collection vials.
- Pour sample into FAST cartridge and elute sample directly into auto-sampler vials.
- Cap vials and put directly onto LC/MS for analysis.

II. FAST Method – Benzodiazepines and Basic Compounds

Sample Dilution Ratio	Sample* Volume	Dilution** Volume
1:1	500 µL	500 µL
1:4	200 µL	800 µL
1:9	100 µL	900 µL





* If sample is hydrolyzed add appropriate aliquot volume after hydrolysis is complete.
** Diluent is 50:50 (Acetonitrile: Distilled Water)

- Sample and diluents are added in an appropriately labeled tube. Add appropriate volume internal standard(s). It is recommended to use an internal standard volume of no more than 200 µL.
- Set up extraction manifold with FAST cartridges and auto-sampler collection vials.
- Pour sample into FAST cartridge and elute sample directly into auto-sampler vials.
- Cap vials and put directly onto LC/MS for analysis.




CONCLUSION

The FAST method outlined is a novel approach to improved sample preparation for LC/MS analysis. The method outlines a simple procedure to prepare urine samples for analysis of multiple drugs and metabolites, by quickly and efficiently reducing the amount of unwanted matrix (through sorbent adsorption) and particulates (filtering through special frits) in the final sample, the analysis can proceed with less chance of matrix suppression and LC column clogging. The FAST method can lengthen the amount of time an LC column can be used for analysis and lower the amount of down time for instrument maintenance. These benefits along with the ability to eliminate the centrifuge and sample transfer steps can lower costs by decreasing turn-around time and reducing instrument and LC column maintenance.

URINE SPECTRA



RESULTS



*NF refers to the 'dilute and shoot' recovery as a normalized referenced (e.g. 100% based on calculated peak areas).
These charts represent 1:1, 1:4 and 1:9 dilution ratios and the % recovery compared to the *NF sample based on calculated peak areas of each compound.

