

**SUMMARY OF SHELF LIFE PRODUCTS**

			Minimum Shelf Life				
		Maximum Shelf Life	guaranteed by		Store at		
Part Number	Description	at Production in months	SKC Ltd in months	Stock held by SKC Ltd	Room Temperature	Store in Freezer	Light Sensitive
225-317	DPM CASSETTE WITH 0.8UM IMPACTOR	12	3	*	*		
225-9001	COATED MCE FILTER & CELLULOSE PAD PRELOADED IN CASSETTE	24	3	*	*		
225-9002	COATED GLASS FILTER &CELLULOSE PAD PRELOADED IN CASSETTE	6	1			*	*
225-9003	COATED GLASS FILTER PRELOADED IN CASSETTE	12	3			*	*
225-9004	COATED GLASS FILTER PRELOADED IN CASSETTE	24	3	*	*		
225-9005	MCE PRE-FILTER&PAD/CELLULOSE FILTER&PAD PRELOADED CASSETTE	12	3	*	*		
225-9006	TEFLON PRE-FILTER & POROUS PLASTIC PAD PRELOADED CASSETTE	12	3	*	*		
225-9007	COATED GLASS FILTER PRELOADED IN CASSETTE	8	3	*		*	
225-9009	COATED GLASS FILTER PRELOADED IN CASSETTE	6	1			*	*
225-9010	COATED GLASS FILTER PRELOADED IN CASSETTE	3	N/A			*	
225-9011	COATED GLASS FILTER TYPE AE	6	1	*		*	
225-9012	COATED GLASS FILTER PRELOADED IN CASSETTE	24	3		*		
225-9013	GLASS FILTER AND CELLULOSE PAD 37MM PRELOADED IN CASSETTE	3	N/A			*	*
225-9014	COATED GLASS FIBRE FILTER X2 PRELOADED IN CASSETTE	1	N/A			*	
225-9015	COATED EXTRA-THICK GLASS FILTER PRELOADED IN CASSETTE	12	3			*	
225-9018	GLASS FIBRE & POLYESTER FILTER 37MM PRELOADED IN CASSETTE	3	N/A			*	
226-10-02	SILICA GEL TUBE 2 PART 10MM X 110MM GS	24	3		*		
226-10-06	SILICA GEL TUBE 2 PART 6MM X 70MM GO	24	3	*	*		
226-104	CHROMOSORB 102 TUBE 2 PART 8MM X 70MM GO	24	3		*		
226-109A	CHROMOSORB 103 TUBE 2 PART 6MM X 85MM GO	24	3		*		
226-117	XAD-2 TUBE 2 PART 6MM X 110MM GS	18	3	*		*	
226-118	XAD-2 TUBE 2 PART 6MM X 110MM GS	18	3	*		*	
226-119	SILICA GEL TUBE 2 PART 6MM X 110MM GS	24	3	*		*	*
226-119A	SILICA GEL TUBE 2 PART 6MM X 110MM GS	24	3	*		*	*
226-120	SILICA GEL TUBE 2 PART 6MM X 110MM GS	24	3			*	*
226-124	PUF/TENAX/PUF TUBE 3 PART 22MM X 100MM GO	24	3		*		
226-125	GLASS BEAD TUBE 7MM X 110MM SILICA GEL TUBE 8MM X 110MM GS	24	3		*		
226-126	PUF/GLASS FIBRE FILTER TUBE 22MM X 100MM GO	24	3		*		
226-127	PORAPAK N/ANASORB CSC CHARCOAL TUBE 2 PART 6MM X 254MM GO	24	3		*		
226-129	PUF/XAD-2/PUF TUBE 3 PART 65MM X 125MM GO	24	3		*		
226-131	PUF TUBE 65MM X 125 MM GO	24	3		*		
226-133	ANASORB 747 TUBE 16MM X 125MM GO	24	3		*		
226-142	CARBON BEADS/TEFLON FILTER TUBE 2PART 16AE 8MM X 85MM GO	24	3		*		
226-143	PUF/XAD-2/PUF TUBE 3 PART 22MM X 100MM GO	24	3		*		
226-151	CHARCOAL TUBE 2 PART 6MM X 70MM GS	12	3			*	
226-153	XAD-2 TUBE 2 PART 8MM X 110MM GS	12	3			*	
226-165	SILICA GEL TUBE	1	N/A			*	
226-171	ANASORB 747/TENAX TA TUBE 16MM X 125MM GO	24	3		*		
226-23	XAD-2 2 PART 6MM X 70MM GS	24	3		*		

**SUMMARY OF SHELF LIFE PRODUCTS**

		Maximum	Shelf Life				
		Shelf Life	guaranteed by		Store at		
		at Production	SKC Ltd	Stock held	Room	Store in	Light
Part Number	Description	in months	in months	by SKC Ltd	Temperature	Freezer	Sensitive
226-27	XAD-2 TUBE 2 PART 8MM X 110MM GS	18	3	*		*	
226-29	ANASORB 747 TUBE 2 PART 8MM X 110MM GS	24	3	*	*		
226-30-07	XAD-2 TUBE 2 PART 8MM X 110MM GS	12	3	*		*	
226-30-16	XAD-2/GLASS FIBRE FILTER TUBE 2PART 13AE 6MM X 75MM GO	24	3	*	*		
226-30-16A	XAD-2/GLASS FIBRE FILTER TUBE 2PART 13AE 6MM X 75MM GO	24	3	*	*		
226-30-18	XAD-2 TUBE 6MM X 70MM GS	24	3	*	*		
226-302	ANASORB CMS TUBE 1/4" X 3" SS	24	3		*		
226-323	TENAX TA TUBE 1/4" X 3" SS	24	3		*		
226-328	ANASORB GCB1 TUBE 1/4" X 3" SS	24	3		*		
226-330	ANASORB CMS/GCB1/GCB2 TD TUBE 3 PART 6MM X 115MM GO	24	3		*		
226-332	TENAX GR TUBE 1/4" X 3" SS	24	3		*		
226-336	TENAX TA TUBE 1/4" X 7" SS	24	3		*		
226-337	ANASORB CMS TUBE 1/4" X 7" SS	24	3		*		
226-338	ANASORB CMS/GCB1/GCB2 TUBE 3 PART 1/4" X 7" SS	24	3		*		
226-339	TENAX TA TD TUBE 14" X 3.5" GO	24	3		*		
226-340	TENAX TA TUBE 1/4" X 3.5" SS	24	3		*		
226-341	ANASORB CMS TUBE 1/4" X 3.5" SS	24	3		*		
226-345	TENAX GR/ANASORB GCB1 TUBE 2 PART 1/4" X 3.5" GO	24	3		*		
226-346	ANASORB GCB1/CMS TUBE 2 PART 1/4" X 3.5" GO	24	3		*		
226-347	ANASORB GCB2/GCB1/CMS TUBE 3 PART 1/4" X 3.5" GO	24	3		*		
226-348	TENAX GR/ANASORB GCB1 TD TUBE 2 PART 1/4" X 3.5" SS	24	3		*		
226-349	ANASORB GCB1/CMS TUBE 2 PART 1/4" X 3.5" SS	24	3		*		
226-350	ANASORB GCB2/GCB1/CMS TUBE 3 PART 1/4" X 3.5" SS	24	3		*		
226-35-01	TENAX TUBE 2 PART 6MM X 70MM GO	24	3		*		
226-35-02	TENAX TUBE 6MM X 130MM GO	24	3		*		
226-352	ANASORB CMS/GCB1/TENAX GR TUBE 1/4" X 4.5" GO	24	3		*		
226-356	ANASORB GCB1 TUBE 1/4" X 3.5" SS	24	3		*		
226-357	TENAX TA TUBE 1/4" X 3.5" SS	24	3		*		
226-358	CHROMOSORB 106 TUBE 2 PART 1/4" X 3.5" SS	24	3		*		
226-359	ANASORB GCB1 TUBE 1/4" X 3.5" GO	24	3		*		
226-360	TENAX TA TUBE 1/4" X 3.5" GO	24	3		*		
226-361	CHROMOSORB 106 TUBE 2 PART 1/4" X 3.5" GO	24	3		*		
226-38-03	CHARCOAL PETROLEUM TUBE 2 PART 6MM X 70MM GS	12	3	*		*	
226-40	MOLECULAR SIEVE TUBE 7X70MM(2) GS & OXIDIZER TUBE 7 X 110MM	24	3	*		*	
226-40-02	MOLECULAR SIEVE TUBE 2 PART 7MM X 110MM GS	24	3	*	*		
226-42	SILICA GEL TUBE 2 PART 8MM X 110MM GS	24	3		*		
226-42-02	FIREBRICK TUBE 7MM X 70MM GS	24	3		*		
226-53	SILICA GEL TUBE 2 PART 6MM X 50MM GO	24	3		*		
226-54	XAD-2 TUBE 2 PART 6MM X 70MM GS	18	3			*	
226-55	SILICA GEL TUBE 2 PART 7MM X 70MM GS	24	3		*		
226-56	TENAX/GLASS FIBRE FILTER TUBE 2PART 13AE 6MM X 75MM GO	24	3	*	*		
226-57	XAD-7/GLASS FIBRE FILTER TUBE 2PART 13AE 6MM X 75MM GO	24	3	*	*		

**SUMMARY OF SHELF LIFE PRODUCTS**

Part Number	Description	Maximum	Shelf Life	Stock held by SKC Ltd	Store at	Store in Freezer	Light Sensitive
		Shelf Life at Production in months	guaranteed by SKC Ltd in months		Room Temperature		
226-57A	XAD-7/GLASS FIBRE FILTER TUBE 2PART 13AE 6MM X 75MM GO	24	3		*		
226-58	XAD-2/QUARTZ FILTER TUBE 2PART 13AE 6MM X 75MM GO	24	3	*	*		
226-58A	XAD-2/QUARTZ FILTER TUBE 2PART 13AE 6MM X 75MM GO	24	3		*		
226-61	SILICA GEL/CHARCOAL TUBE 3 PART 10MM X 210MM GS	24	3		*		
226-63	SILICA GEL TUBE 6MM X 40MM GO	24	3		*		
226-67	ANASORB CSC COCONUT CHARCOAL TUBE 2 PART 6MM X 70MM GS	24	3		*		
226-68	JXC CHARCOAL/DRIERITE TUBE 3PART 8AE 6MM X 160MM GS	12	3			*	
226-70A	SILICA GEL TUBE 2 PART 8MM X 150MM GS	12	3	*		*	
226-73	ANASORB CSC COCONUT CHARCOAL TUBE 2 PART 6MM X 70MM GS	24	3	*	*		
226-80	ANASORB 747 TUBE 2 PART 6MM X 70MM GS	24	3	*	*		
226-92	PUF TUBE 22MM X 100MM GO	24	3	*	*		
226-96	XAD-7 TUBE 2 PART 8MM X 110MM GS	24	3	*	*		
226-98	XAD-7 TUBE 2 PART 6MM X 70MM GS.	24	3	*	*		
226-99	SILICA GEL/G F FILTER OVS TUBE 2 PART 13\6MM X 75MM GO	24	3	*	*		
227-101A	ANACHEK PRESPIKED TUBE 2 PART FOR TOLUENE	12	3			*	
227-101B	ANACHEK PRESPIKED TUBE 2 PART FOR TOLUENE	12	3			*	
227-102A	ANACHEK PRESPIKED TUBE 2 PART FOR TRICHLOROETHYLENE	12	3			*	
227-102B	ANACHEK PRESPIKED TUBE 2 PART FOR TRICHLOROETHYLENE	12	3			*	
227-103A	ANACHEK PRESPIKED TUBE 2 PART FOR ETHYL ACETATE	12	3			*	
227-103B	ANACHEK PRESPIKED TUBE 2 PART FOR ETHYL ACETATE	12	3			*	
227-104A	ANACHEK PRESPIKED TUBE 2 PART FOR METHYLENE CHLORIDE	12	3			*	
227-104B	ANACHEK PRESPIKED TUBE 2 PART FOR METHYLENE CHLORIDE	12	3			*	
227-105A	ANACHEK PRESPIKED TUBE 2 PART FOR HEXANE	12	3			*	
227-105B	ANACHEK PRESPIKED TUBE 2 PART FOR HEXANE	12	3			*	
227-106A	ANACHEK PRESPIKED TUBE 2 PART FOR BUTYL ACETATE	12	3			*	
227-106B	ANACHEK PRESPIKED TUBE 2 PART FOR BUTYL ACETATE	12	3			*	
227-501	ANACHEK PRESPIKED PASSIVE SAMPLER FOR TOLUENE	6	1			*	
227-502	ANACHEK PRESPIKED PASSIVE SAMPLER FOR TRICHLOROETHYLENE	6	1			*	
227-503	ANACHEK PRESPIKED PASSIVE SAMPLER FOR ETHYL ACETATE	6	1			*	
227-504	ANACHEK PRESPIKED PASSIVE SAMPLER FOR METHYLENE CHLORIDE	3	N/A			*	
227-505	ANACHEK PRESPIKED PASSIVE SAMPLER FOR HEXANE	6	1			*	
227-506	ANACHEK PRESPIKED PASSIVE SAMPLER FOR BUTYL ACETATE	6	1			*	
500-100	UMEX100 PASSIVE SAMPLER FOR FORMALDEHYDE	6	1	*		*	*
520-02A	SORBENT CAPSULE	24	3	*	*		
530-11	ANASORB CA CAPSULES	24	3	*	*		
530-14	ANASORB LS CAPSULES	24	3		*		
575-001	PASSIVE SAMPLER FOR ORGANIC VAPOUR	24	3	*	*		
575-001A	PASSIVE SAMPLER FOR ORGANIC VAPOUR	24	3	*	*		
575-001MC	PASSIVE SAMPLER FOR METHYLENE CHLORIDE	24	3		*		
575-001MCA	PASSIVE SAMPLER FOR METHYLENE CHLORIDE	24	3		*		

**SUMMARY OF SHELF LIFE PRODUCTS**

		<b>Maximum Shelf Life at Production</b>	<b>Shelf Life guaranteed by SKC Ltd</b>	<b>Stock held by SKC Ltd</b>	<b>Store at Room Temperature</b>	<b>Store in Freezer</b>	<b>Light Sensitive</b>
<b>Part Number</b>	<b>Description</b>	<b>in months</b>	<b>in months</b>				
575-002	PASSIVE SAMPLER FOR ORGANIC VAPOUR	24	3	*	*		
575-003	PASSIVE SAMPLER FOR ORGANIC VAPOUR	24	3	*	*		
575-005	PASSIVE SAMPLER FOR ETHYLENE OXIDE	12	3			*	
P226129R	REPACK OF 226-129	24	3		*		
P226131R	REPACK OF 226-131	24	3		*		
P226143R	REPACK OF 226-143	24	3		*		
P22692R	REPACK OF 226-92	24	3		*		